Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-COLLAGEN TYPE X A1 (ZEBRAFISH) ANTIBODY

SKU: EB07975

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

collagen type X alpha 1|wu:fb10c08|wu:fa66d11|fb10c08|fa66d11|col10a1|collagen, type X, alpha 1 Alias

Names

Accession NP 001077296.1 ID

Blocking

EBP07975 **Peptide**

Peptide with sequence C-YNAEQHYDPET, from the internal region of the protein sequence according to Immunogen

NP 001077296.1.

Peptide

C-YNAEQHYDPET Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Zebrafish

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection Antibody detection limit dilution 1:16000.

Limit

Western Blot: Not yet tested. At this stage we are dependent on researchers in the field for further

Western

Blot

characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported

with other antibodies/lysates?

Application

Pep-ELISA **Type**