


GOAT ANTI-CONNECTIVE TISSUE GROWTH FACTOR ANTIBODY

SKU: EB06683



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias connective tissue growth factor|IGFBP8|NOV2|CCN2|CTGF

Names

Usage Summary **Immunofluorescence:** Strong expression of the protein seen in the cytoplasm of NIH3T3 cells. Recommended concentration: 5µg/ml.

Accession ID NP_001892

Blocking Peptide EBP06683

Immunogen Peptide with sequence C-SLYYRKMYGDMA, from the C Terminus of the protein sequence according to NP_001892.

Peptide Sequence C-SLYYRKMYGDMA

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Mouse, Rat
Human Gene ID	1490
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Rat Liver. Recommended concentration: 6-8µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:1000.
Application Type	Pep-ELISA, IF, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

