

GOAT ANTI-CORO2B / CLIPINC ANTIBODY

SKU: EB06362



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	coronin, actin-binding protein, 2B coronin, actin-binding, 2B clipin C coronin, actin binding protein, 2B KIAA0925 ClipinC CLIPINC CORO2B
Accession ID	NP_006082.3 ; NP_001177385.1
Blocking Peptide	EBP06362
Immunogen	Peptide with sequence C-QLQLELKNLRNSPK, from the C Terminus of the protein sequence according to NP_006082.3 ; NP_001177385.1.
Product Comments	This epitope was chosen to prevent cross-reactivity with hypothetical protein XP_121721.
Peptide Sequence	C-QLQLELKNLRNSPK
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	10391
Mouse Gene ID	235431
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Preliminary experiments gave no signal but low background in human brain and testis lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)