Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-CORO2B / CLIPINC ANTIBODY

SKU: EB06362

SPECIFICATIONS

Unit Size

 $100 \mu g$

Synonym / Alias

coronin, actin-binding protein, 2B|coronin, actin-binding, 2B|clipin C|coronin, actin binding protein,

Names

2B|KIAA0925|ClipinC|CLIPINC|CORO2B

Accession

ID

NP_006082.3; NP_001177385.1

Blocking

Peptide

EBP06362

Immunogen

Peptide with sequence C-QLQLELKNLRNSPK, from the C Terminus of the protein sequence according to

NP_006082.3; NP_001177385.1.

Product

Comments

This epitope was chosen to prevent cross-reactivity with hypothetical protein XP_121721.

Peptide

Sequence

C-QLQLELKNLRNSPK

Shipping

Instructions Refrigerated

Predicted

Species

Human, Mouse

Human

Gene ID

10391

Mouse

235431

Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection Limit

Antibody detection limit dilution 1:16000.

Western

Blot

Preliminary experiments gave no signal but low background in human brain and testis lysates at up to 1µg/ml. We

would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?

Application

Type

Pep-ELISA