

GOAT ANTI-CORONIN 3 / CORONIN 1C ANTIBODY

SKU: EB06357

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	coronin, actin-binding protein, 1C coronin 1C coronin, actin binding protein, 1C HCRNN4 CORO1C coronin 1C coronin-3
Accession ID	NP_055140.1
Blocking Peptide	EBP06357
Immunogen	Peptide with sequence C-RISKLEQQMAKIAA, from the C Terminus of the protein sequence according to NP_055140.1.
Peptide Sequence	C-RISKLEQQMAKIAA
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Pig, Cow
Reactive Species	Human
Human Gene ID	23603
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 55kDa band observed in human lung and placenta lysates (calculated MW of 53.2kDa according to NP_055140.1). Recommended for use at 0.2-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

