

GOAT ANTI-COUP-TF I / EAR3 / NR2F1 ANTIBODY

SKU: EB05637



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	transcription factor COUP 1 (chicken ovalbumin upstream promoter 1, v-erb-a homolog-like 3) Transcription factor COUP 1 (chicken ovalbumin upstream promoter 1, TCFCOUP1 COUP-
Names	TFI TFCOUP1 ERBAL3 SVP44 NR2F2 EAR-3 EAR3 nuclear receptor subfamily 2, group F, member 1 NR2F1
Accession ID	NP_005645.1
Blocking Peptide	EBP05637
Immunogen	Peptide with sequence AMVVSSWRDPQDD-C, from the N Terminus of the protein sequence according to NP_005645.1.
Peptide Sequence	AMVVSSWRDPQDD-C
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Human Gene ID	7025
Mouse Gene ID	13865
Rat Gene ID	81808
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	No signal obtained yet but low background observed in HeLa , HepG2 , Human Brain, Human Heart lysates at upto 3µg/ml.
Application Type	Pep-ELISA