

GOAT ANTI-CPT1A (MOUSE) ANTIBODY

SKU: EB08924



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	carnitine palmitoyltransferase 1, liver carnitine palmitoyltransferase 1A C730027G07 L-CPT1 CPT1-L CPT1 carnitine palmitoyltransferase 1A, liver Cpt1a
Accession ID	NP_038523.2
Blocking Peptide	EBP08924
Immunogen	Peptide with sequence C-PETDSHRFGKHLRQ, from the C Terminus of the protein sequence according to NP_038523.2.
Peptide Sequence	C-PETDSHRFGKHLRQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Mouse, Rat, Pig
Mouse Gene ID	12894
Rat Gene ID	25757
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)