

GOAT ANTI-CREB3 (AA76-87) ANTIBODY

SKU: EB12173

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	CREB3 cAMP responsive element binding protein 3 LUMAN LZIP CREB-3 basic leucine zipper protein cAMP-responsive element-binding protein 3 cyclic AMP response element (CRE)-binding protein/activating transcription factor 1 cyclic AMP-responsive element-binding protein 3 leucin zipper protein transcription factor LZIP-alpha
Accession ID	NP_006359.3
Blocking Peptide	EBP12173
Immunogen	Peptide with sequence C-HHDHTYSLPRET, from the internal region of the protein sequence according to NP_006359.3.
Peptide Sequence	C-HHDHTYSLPRET
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	10488
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Skin. Recommended concentration: 5µg/ml.
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 45kDa band observed in lysates of cell line HepG2 (calculated MW of 41.4kDa according to NP_006359.3). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

