Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-KAT3A / CBP ANTIBODY

SKU: EB07833

SPECIFICATIONS

Unit Size $100 \mu g$

Synonym / Alias

KAT3A|CREB-binding protein|CREB binding protein|RTS|RSTS|CBP|CREB binding protein (Rubinstein-Taybi

syndrome)|CREBBP Names

Accession

ID

NP_004371.2; NP_001073315.1

Blocking

Peptide

EBP07833

Peptide with sequence C-DTTGDTLEKFVEG, from the C Terminus of the protein sequence according to **Immunogen**

NP_004371.2; NP_001073315.1.

Product

This antibody is expected to recognise both reported isoforms (NP_004371.2 and NP_001073315.1). **Comments**

Peptide

C-DTTGDTLEKFVEG Sequence

Shipping

Instructions Refrigerated

Predicted

Species

Human, Mouse, Rat, Dog

Human

1387 Gene ID

Mouse Gene ID

12914

Rat Gene ID 54244

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/aspiring_medium.png$

Grade **ELISA**

Detection

Antibody detection limit dilution 1:64000.

Limit

Western Blot

Preliminary experiments in lysates of human breast cancer and of cell lines HeLa and A431 gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in

the field - have any results been reported with other antibodies/lysates?

Application

Type

Pep-ELISA