



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600


## **GOAT ANTI-CRYAB ANTIBODY**

**SKU:** EB08149



Telephone: (650) 697-3600







Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

**Alias** heat-shock 20 kD like-protein|HSPB5|CTPP2|CRYA2|crystallin, alpha B|CRYAB

Names

**Accession** 

NP\_001876.1 ID

**Blocking** 

EBP08149 **Peptide** 

Peptide with sequence C-RLEKDRFSVNLD, from the internal region of the protein sequence according to **Immunogen** 

NP 001876.1.

**Peptide** Sequence

C-RLEKDRFSVNLD

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat **Species** 

Reactive

Mouse, Rat **Species** 

Human

1410 Gene ID

Mouse

12955 **Gene ID** 

Rat Gene ID 25420

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**Grade ELISA** 

**Detection** Antibody detection limit dilution 1:128000.

Limit

Western Approx 22kDa band observed in Mouse and Rat Eye lysates (calculated MW of 20.1kDa according to Mouse

NP\_034094.1 and Rat NP\_037067.1). Recommended concentration: 0.03-0.1-µg/ml. Blot

**Application** 

Pep-ELISA, WB

**Type** 





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GALLERY IMAGES**

