

GOAT ANTI-CST3 / CYSTATIN C ANTIBODY

SKU: EB09287

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	post-gamma-globulin neuroendocrine basic polypeptide gamma-trace cystatin C (amyloid angiopathy and cerebral hemorrhage) cystatin 3 bA218C14.4 (cystatin C) OTTHUMP00000164182 OTTHUMP00000164181 OTTHUMP00000030440 MGC117328 ARMD11 cystatin C CST3
Accession ID	NP_000090.1
Blocking Peptide	EBP09287
Immunogen	Peptide with sequence C-SVEEEGVRRALD, from the internal region of the protein sequence according to NP_000090.1.
Peptide Sequence	C-SVEEEGVRRALD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	1471
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx. 15kDa band observed in Human Testis and Prostate lysates (calculated MW of 15.8kDa according to NP_000090.1). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa