


GOAT ANTI-CXCR6 ANTIBODY

SKU: EB09999



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names CXCR6|TYMSTR|STRL33|OTTHUMP00000209998|OTTHUMP00000209997|OTTHUMP00000209996|G protein-coupled receptor TYMSTR|G protein-coupled receptor|chemokine (C-X-C motif) receptor 6|CD186|BONZO

Accession ID NP_006555.1

Blocking Peptide EBP09999

Immunogen Peptide with sequence C-SEDNSKTFASAHN, from the C Terminus of the protein sequence according to NP_006555.1.

Peptide Sequence C-SEDNSKTFASAHN

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	10663
Mouse Gene ID	80901
Rat Gene ID	100124593
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Thymus. Recommended concentration: 3.75µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx 40-45kDa band observed in HEK293 lysates (calculated MW of 39.3kDa according to NP_006555.1). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

