

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-CXCR4B (ZEBRAFISH) ANTIBODY

SKU: EB08008

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias cxcr4b|chemokine (C-X-C motif), receptor 4b|CXCR4|cb403|zgc:109863|ody|odysseus

Names

Accession NP 571909.1 ID

Blocking EBP08008

Peptide

Immunogen Peptide with sequence CELTYPQEAN, from the internal region of the protein sequence according to NP 571909.1.

Product

This antibody is not expected to cross-react with Cxcr4a. Comments

Peptide

CELTYPOEAN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Zebrafish

Species Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

FLISA

Detection Antibody detection limit dilution 1:8000.

Limit

Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We

Blot would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?

Application

Western

Pep-ELISA **Type**