

GOAT ANTI-CYB561D2 ANTIBODY

SKU: EB05655



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	OTTHUMP0000021067 OTTHUMP00000210676 OTTHUMP00000210675 TSP10 putative tumor suppressor 101F6 cytochrome b-561 domain containing 2 CYB561D2
Accession ID	NP_008953.1
Blocking Peptide	EBP05655
Immunogen	Peptide with sequence C-VSNAYLYRKRIQP, from the C Terminus of the protein sequence according to NP_008953.1.
Peptide Sequence	C-VSNAYLYRKRIQP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	11068
Mouse Gene ID	56368
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 25kDa band observed in human liver lysates (predicted MW of 26kDa according to NP_008953). Recommended for use at 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

