

GOAT ANTI-CYCLIN D1 ANTIBODY

SKU: EB12329

250kDa

150kDa

100kDa

75kDa


50kDa

37kDa

25kDa

20kDa

15kDa



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PRAD1 oncogene G1/S-specific cyclin-D1 BCL-1 oncogene B-cell lymphoma 1 protein B-cell CLL/lymphoma 1 U21B31 PRAD1 D11S287E BCL1 cyclin D1 CCND1
Accession ID	NP_444284.1
Blocking Peptide	EBP12329
Immunogen	Peptide with sequence TRFLSRVIKCDPD, from the internal region of the protein sequence according to NP_444284.1.
Peptide Sequence	TRFLSRVIKCDPD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	595
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 30kDa band observed in Human Placenta lysates (calculated MW of 33.7kDa according to NP_444284.1). Recommended concentration: 0.5-1.5µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

