

GOAT ANTI-CYP2B6 (AA247-259) ANTIBODY

SKU: EB12714



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	cytochrome P450, subfamily IIB (phenobarbital-inducible), polypeptide 6 cytochrome P450 IIB1 cytochrome P450 2B6 1, 4-cineole 2-exo-monooxygenase P450 IIB1 EFVM CYP11B6 CYP2B7P CYP2B7 CYP2B CPB6 cytochrome P450, family 2, subfamily B, polypeptide 6 CYP2B6
Accession ID	NP_000758.1
Blocking Peptide	EBP12714
Immunogen	Peptide with sequence C-HSVEKHRETLDPS, from the internal region of the protein sequence according to NP_000758.1.
Peptide Sequence	C-HSVEKHRETLDPS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	1555
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 75kDa band observed in lysates of cell line HepG2 and K562 (calculated MW of 56.3kDa according to NP_000758.1). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

