


GOAT ANTI-CYP2D5 / CYP2D10 (RAT) ANTIBODY

SKU: EB10750



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	cytochrome P450CMF cytochrome P450-DB5 cytochrome P450, family 2, subfamily d, polypeptide 10 cytochrome P450 2D10 P450-DB5 P450-CMF1B CYP11D10 MGC108540 Cyp2d10 cytochrome P450, family 2, subfamily d, polypeptide 5 Cyp2d5
Accession ID	NP_775426.1
Blocking Peptide	EBP10750
Immunogen	Peptide with sequence CDAFTAQNGRSINPK, from the internal region of the protein sequence according to NP_775426.1.
Peptide Sequence	CDAFTAQNGRSINPK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Rat
Reactive Species	Rat
Rat Gene ID	286963
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 65kDa band observed in Rat Spleen lysates (calculated MW of 57.1kDa according to NP_775426.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

