

GOAT ANTI-CYSTATIN F / CMAP ANTIBODY

SKU: EB05429



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	dJ568C11.1 (cystatin F (leukocystatin)) cystatin F cystatin-like metastasis-associated protein leukocystatin cystatin 7 CMAP cystatin F (leukocystatin) CST7
Accession ID	NP_003641.2
Blocking Peptide	EBP05429
Immunogen	Peptide with sequence PWLQHFEVPVLRCH, from the C Terminus of the protein sequence according to NP_003641.2.
Peptide Sequence	PWLQHFEVPVLRCH
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Mouse
Human Gene ID	8530
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 22kDa band seen in NSO lysate. Recommended for use at 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

