

GOAT ANTI-DAB2 / DOC2 (AA 345-358) ANTIBODY

SKU: EB10685



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	mitogen-responsive phosphoprotein disabled homolog 2 differentially-expressed protein 2 FLJ26626 DOC2 DOC-2 disabled homolog 2, mitogen-responsive phosphoprotein (Drosophila) DAB2
Accession ID	NP_001334.2
Blocking Peptide	EBP10685
Immunogen	Peptide with sequence C-QQFDQISNRTGKQE, from the internal region of the protein sequence according to NP_001334.2.
Peptide Sequence	C-QQFDQISNRTGKQE
Shipping Instructions	Refrigerated
Predicted Species	Human, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	1601
Rat Gene ID	79128
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 100kDa band observed in lysates of cell line HeLa (calculated MW of 82.4kDa according to NP_001334.2). Recommended concentration: 0.03-0.1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

