

GOAT ANTI-DAB2 / DOC2 ANTIBODY

SKU: EB06077

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	FLJ26626 mitogen-responsive phosphoprotein disabled homolog 2 disabled (Drosophila) homolog 2 (mitogen-responsive phosphoprotein) mitogen-responsive phosphoprotein DOC-2 disabled homolog 2, mitogen-responsive phosphoprotein (Drosophila) DOC2 DAB2
Accession ID	NP_001334.2
Blocking Peptide	EBP06077
Immunogen	Peptide with sequence C-SEMYRDPFGNPFA, from the C Terminus of the protein sequence according to NP_001334.2.
Peptide Sequence	C-SEMYRDPFGNPFA
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	1601
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 75kDa band observed in human kidney lysates (calculated MW of 82.6kDa according to NP_001334.1). Recommended concentration: 0.01-0.03µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

