

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-DAB2 / DOC2 ANTIBODY

SKU: EB06077







250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Unit Size 100 µg

Synonym / FLJ26626|mitogen-responsive phosphoprotein|disabled homolog 2|disabled (Drosophila) homolog 2 (mitogen-

Alias responsive phosphoprotein)|mitogen-responsive phosphoprotein|DOC-2|disabled homolog 2, mitogen-

responsive phosphoprotein (Drosophila)|DOC2|DAB2 **Names**

Accession

NP_001334.2

Blocking

EBP06077 **Peptide**

Peptide with sequence C-SEMYRDPFGNPFA, from the C Terminus of the protein sequence according to **Immunogen**

NP 001334.2.

Peptide

C-SEMYRDPFGNPFA Sequence

Refrigerated Instructions

Predicted

Species

Human

Reactive

Human **Species**

Human

1601 Gene ID

Product

Grade ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western Approx 75kDa band observed in human kidney lysates (calculated MW of 82.6kDa according to NP_001334.1).

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png

Blot Recommended concentration: 0.01-0.03µg/ml.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa