

GOAT ANTI-DAPK1 ANTIBODY

SKU: EB05274



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	DAPK death-associated protein kinase 1 DAPK1
Accession ID	NP_004929.1 GI:4826684
Immunogen	Peptide with sequence C-SGTSYNSISSVVSRR, from the C Terminus of the protein sequence according to NP_004929.1 GI:4826684.
Peptide Sequence	C-SGTSYNSISSVVSRR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	1612
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Mouse and Human Brain extracts and HeLa cell lysates at upto 4µg/ml.