



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-DAPK1 ANTIBODY

SKU: EB05274

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias DAPK|death-associated protein kinase 1|DAPK1

Names

Accession NP 004929.1 GI:4826684

Peptide with sequence C-SGTSYNSISSVVSR, from the C Terminus of the protein sequence according to **Immunogen** NP_004929.1 GI:4826684.

Peptide

C-SGTSYNSISSVVSR Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

1612 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

Western Western blot: No signal obtained yet but low background observed in Mouse and Human Brain extracts and HeLa

Blot cell lysates at upto 4µg/ml.