Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-DBF4B ANTIBODY

SKU: EB10066

SPECIFICATIONS

Unit Size 100 µg

Synonym / DBF4B|zinc finger, DBF-type containing 1B|ZDBF1B|MGC15009|FLJ13087|DRF1|Dbf4-related factor 1|DBF4 Alias homolog B (S. cerevisiae)|DBF4 homolog B|chiffon homolog B|chifb|ASKL1|activator of S-phase kinase-like protein

Names

Accession

ID

NP_663696.1; NP_079380.1

Blocking EBP10066

Peptide

Peptide with sequence C-PHDTTPLHEEVSP, from the internal region of the protein sequence according to **Immunogen**

NP 663696.1; NP 079380.1.

Product

This antibody is expected to recognize reported isoform 1 (NP_663696.1) only. Comments

Peptide

C-PHDTTPLHEEVSP Sequence

Refrigerated Instructions

Predicted

Human **Species**

Human 80174

Gene ID

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:4000.

Limit

Western Blot

Preliminary experiments in Human Testis lysates gave no specific signal but low background (at antibody concentration up to $1\mu g/ml$). We would appreciate any feedback from people in the field - have any results been

reported with other antibodies/lysates?

Application

Pep-ELISA **Type**