

GOAT ANTI-DBF4B ANTIBODY

SKU: EB10066



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	DBF4B zinc finger, DBF-type containing 1B ZDBF1B MGC15009 FLJ13087 DRF1 Dbf4-related factor 1 DBF4 homolog B (S. cerevisiae) DBF4 homolog B chiffon homolog B chifb ASKL1 activator of S-phase kinase-like protein 1
Accession ID	NP_663696.1; NP_079380.1
Blocking Peptide	EBP10066
Immunogen	Peptide with sequence C-PHDTTPLHEEVSP, from the internal region of the protein sequence according to NP_663696.1; NP_079380.1.
Product Comments	This antibody is expected to recognize reported isoform 1 (NP_663696.1) only.
Peptide Sequence	C-PHDTTPLHEEVSP
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	80174
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Preliminary experiments in Human Testis lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

