


GOAT ANTI-DDX58 ANTIBODY

SKU: EB08003



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	retinoic acid-inducible gene 1 RNA helicase RIG-1 RNA helicase OTTHUMP00000045225 DEAD/H (Asp-Glu-Ala-Asp/His) box polypeptide RIG-1 DEAD/H (Asp-Glu-Ala-Asp/His) box polypeptide RIG-1 FLJ13599 DKFZp686N19181 DKFZp434J1111 DEAD (Asp-Glu-Ala-Asp) box polypeptide 58 DDX58
Usage Summary	Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 and HeLa cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of U937 cells. Recommended concentration: 10ug/ml.
Accession ID	NP_055129.2
Blocking Peptide	EBP08003
Immunogen	Peptide with sequence C-KDFHFEKIPFDPAE, from the C Terminus of the protein sequence according to NP_055129.2.
Peptide Sequence	C-KDFHFEKIPFDPAE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Human
Human Gene ID	23586
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png
IHC Results	Paraffin embedded Human Spleen, Placenta, and Brain (Cortex). Recommended concentration: 3.75µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, FC, IF, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

