

EBP06822 - IMMUNIZING PEPTIDE

SKU: EBP06822



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_922938.1; NP_922939.1 NP_922943
Immunogen	Peptide with sequence C-SRNKENHSQPTQ, from the internal region of the protein sequence according to NP_922938.1; NP_922939.1 NP_922943.
Product Comments	The peptide was used in the production of EBP06822 - Goat Anti-DECTIN1 / beta-Glucan receptor Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-SRNKENHSQPTQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Preliminary experiments in lysates of human liver, bone marrow and peripheral blood mononucleocytes gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

DOCUMENTS

- [Data Sheet](#)