

GOAT ANTI-DECTIN1 / BETA-GLUCAN RECEPTOR ANTIBODY

SKU: EB06822



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	dendritic cell-associated C-type lectin 1 dectin-1 transmembrane protein dectin-1 lectin-like receptor 1 dendritic cell-associated C-type lectin-1 beta-glucan receptor C-type lectin, superfamily member 12 C-type (calcium dependent, carbohydrate-recognition domain) lectin, superfamily member 12 DECTIN1 CLECSF12 BGR HGNC:14558 C-type lectin domain family 7, member A CLEC7A
Accession ID	NP_922938.1; NP_922939.1 NP_922943
Blocking Peptide	EBP06822
Immunogen	Peptide with sequence C-SRNKENHSQPTQ, from the internal region of the protein sequence according to NP_922938.1; NP_922939.1 NP_922943.
Product Comments	This antibody is expected to recognise isf a (NP_922938), isf c (NP_922939), and isf g (NP_922943)
Peptide Sequence	C-SRNKENHSQPTQ
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog
Human Gene ID	64581
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Preliminary experiments in lysates of human liver, bone marrow and peripheral blood mononucleocytes gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA