
Telephone: (650) 697-3600





## **GOAT ANTI-DEFB107A / DEFB107B ANTIBODY**

**SKU:** EB11110



## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

BD-7| beta-defensin 107| beta-defensin 7| DEFB107| DEFB107B| DEFB-7| defensin, beta 107A| defensin, beta

7|DEFB107A **Names** 

Accession

Alias

NP\_001032757.2; NP\_001035795.1 ID

**Blocking** 

EBP11110 **Peptide** 

Peptide with sequence HRALISKRMEGHCE, from the internal region of the protein sequence according to **Immunogen** 

NP\_001032757.2; NP\_001035795.1.

**Peptide Sequence** 

HRALISKRMEGHCE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

**Predicted** 

Human

**Species** Reactive

Human

**Species** 

Human **Gene ID** 

245910, 503614

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

IHC Results Paraffin embedded Human Colon. Recommended concentration: 5µg/ml.

**ELISA** 

**Detection** Antibody detection limit dilution 1:32000.

Limit

Blot

Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as

the calculated size of 7.85kDa according to NP\_001032757.2.

Application

Western

Pep-ELISA, IHC Type

## GALLERY IMAGES

