



GOAT ANTI-DEFENSIN, BETA 1 (MOUSE) ANTIBODY

SKU: EB11073

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Defb1|beta-defensin 1|AW260221| defensin beta 1| BD-1| AL590630.1

Names

Accession NP 031869.1 ID

Blocking

EBP11073 **Peptide**

Immunogen Peptide with sequence TKLQGTCKPDKPN, from the C Terminus of the protein sequence according to NP 031869.1.

Product

ADJUST NAMES AND SYMBOLS Comments

Peptide

TKLOGTCKPDKPN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Instructions Refrigerated

Predicted

Mouse, Rat **Species**

Mouse

13214 **Gene ID** Rat Gene ID 83687

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:4000.

Limit

Application Pep-ELISA

Type