

GOAT ANTI-DEFENSIN, BETA 1 (MOUSE) ANTIBODY

SKU: EB11073



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Defb1 beta-defensin 1 AW260221 defensin beta 1 BD-1 AL590630.1
Accession ID	NP_031869.1
Blocking Peptide	EBP11073
Immunogen	Peptide with sequence TKLQGTCKPDKPN, from the C Terminus of the protein sequence according to NP_031869.1.
Product Comments	ADJUST NAMES AND SYMBOLS
Peptide Sequence	TKLQGTCKPDKPN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Mouse, Rat
Mouse Gene ID	13214
Rat Gene ID	83687
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)