

GOAT ANTI-DIAPH3 ANTIBODY

SKU: EB08152

SPECIFICATIONS

Unit Size	100 µg
Synonym /	diaphanous-related formin 3 diaphanous homolog
Alias	3 diap3 FLJ34705 DRF3 DKFZp686A13178 DKFZP434C0931 RP11-26P21.1 diaphanous homolog 3
Names	(Drosophila) DIAPH3
Accession ID	NP_001035982.1; NP_112194.2; NP_001245295.1; NP_001245296.1; NP_001245297.1; NP_001245298.1
Blocking Peptide	EBP08152
Immunogen	Peptide with sequence C-RTPMPKDVRQS, from the internal region (near the C terminus) of the protein sequence according to NP_001035982.1; NP_112194.2; NP_001245295.1; NP_001245296.1; NP_001245297.1; NP_001245298.1.
Peptide Sequence	C-RTPMPKDVRQS
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	81624
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)