


GOAT ANTI-DISHEVELLED1 (AA20-32) ANTIBODY

SKU: EB11466



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names segment polarity protein dishevelled homolog DVL-1| dishevelled-1| DSH homolog 1| mKIAA4029| Dvl| RP23-242K15.15-003| dishevelled, dsh homolog 1 (Drosophila)|Dvl1

Accession ID NP_034221.3

Blocking Peptide EBP11466

Immunogen Peptide with sequence KLPVAPERVTAD-C, from the N Terminus of the protein sequence according to NP_034221.3.

Peptide Sequence KLPVAPERVTAD-C

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human, Mouse, Rat, Pig, Cow
Reactive Species	Human
Human Gene ID	1855
Mouse Gene ID	13542
Rat Gene ID	83721
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Skeletal Muscle and Kidney. Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA, IHC

SELECTED REFERENCES

[{"pmid": 24739947, "intro": "", "title": "Dickkopf-related protein 1 inhibits the WNT signaling pathway and improves pig\r\nnoocyte maturation.", "author": "Spate LD, Brown AN, Redel BK, Whitworth KM, Murphy CN, Prather RS.", "journal": "PLoS One. 2014 Apr 16;9(4):e95114."}]

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

