


GOAT ANTI-DKC1 ANTIBODY

SKU: EB07703



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias XAP101|NOLA4|NAP57|DKC|dyskeratosis congenita 1, dyskerin|DKC1

Names

Accession ID NP_001354.1

Blocking Peptide EBP07703

Immunogen Peptide with sequence C-KRKRESESESDETPP, from the internal region of the protein sequence according to NP_001354.1.

Peptide Sequence C-KRKRESESESDETPP

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human, Mouse
Reactive Species	Human
Human Gene ID	1736
Mouse Gene ID	245474
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Spleen shows strong nuclear staining of select splenocytes. Recommended concentration, 5-10µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

