

GOAT ANTI-DKK1 ANTIBODY

SKU: EB05358



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	dickkopf (Xenopus laevis) homolog 1 dickkopf-1 like DKK-1 SK dickkopf homolog 1 (Xenopus laevis) DKK1
Names	
Accession ID	NP_036374.1
Blocking Peptide	EBP05358
Immunogen	Peptide with sequence HQASNSSLHTCQRH, from the C Terminus of the protein sequence according to NP_036374.1.
Peptide Sequence	HQASNSSLHTCQRH
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Mouse
Human Gene ID	22943
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 35kDa band seen in Mouse Liver lysate. Recommended for use at 1-3µg/ml.