



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-DKK2 ANTIBODY

SKU: EB05360



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa 37kDa

25kDa 20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

dickkopf related protein-2|dickkopf homolog 2|dickkopf 2|dickkopf (Xenopus laevis) homolog 2|Dickkopf gene

2|DKK-2|dickkopf homolog 2 (Xenopus laevis)|DKK2

Names **Accession**

NP_055236.1

Blocking

EBP05360 **Peptide**

Peptide with sequence ATYSSKARLHVCQKI, from the C Terminus of the protein sequence according to

NP 055236.1.

Product From sequence analysis it is possible that this product will crossreact with DKK4, unfortunately this has not yet

been specifically tested. Comments

Peptide

Sequence

Immunogen

ATYSSKARLHVCQKI

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Human

27123

Predicted

Species

Human, Mouse, Rat, Dog, Cow

Reactive

Species

Human

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Detection

Antibody detection limit dilution 1:8000.

Limit

Western Approx 28kDa band observed in Human Heart lysate (calculated MW of 28.5kDa according to NP_055236).

Recommended concentration: 1-3µg/ml Blot

Application

Type

Pep-ELISA, WB

GALLERY IMAGES





250102

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

150kDa
100kDa
75kDa
50kDa
37kDa
25kDa

15kDa

20kDa