



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-DKK4 ANTIBODY

SKU: EB05362



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names DKK4|dickkopf homolog 4 (Xenopus laevis)|DKK-4|hdkk-4|Dickkopf gene 4|dickkopf (Xenopus laevis) homolog

4|MGC129562|MGC129563|dickkopf homolog 4

Accession

ID

NP_055235.1

Blocking

Peptide

EBP05362

Peptide with sequence RQHARLRVCQKIEKL, from the C Terminus of the protein sequence according to **Immunogen**

NP 055235.1.

Peptide

Sequence

ROHARLRVCOKIEKL

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Human

27121

Predicted

Species

Human, Mouse, Rat, Dog, Pig

Reactive

Species

Human

Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection

Antibody detection limit dilution 1:32000.

Limit

Western Blot

Approx 30-35kDa band seen in Jurkat lysate. Recommended for use at 1-3µg/ml. Please note that literature indicates DKK4 is detected as three major immunoreactive species of approx. 40 kDa [form (i)], 30-32 kDa [form (ii)] and 15-17 kDa [form (iii)]. We believe this antibody will recognise 30-32kDa species of DKK4.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa