

GOAT ANTI-DMBX1 / OTX3 / MBX ANTIBODY

SKU: EB05944



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	orthodenticle homolog 3 DMBX1 paired-like homeobox protein DMBX1-like homeoprotein MBX orthodenticle homolog 3 (Drosophila) PAXB MBX OTX3 diencephalon/mesencephalon homeobox 1
Accession ID	NP_671725.1; NP_757379.1
Blocking Peptide	EBP05944
Immunogen	Peptide with sequence C-QHAASLGDLTPN, from the C Terminus of the protein sequence according to NP_671725.1; NP_757379.1.
Product Comments	This antibody is expected to recognise both reported human isoforms (as represented by NP_671725.1 and NP_757379.1).
Peptide Sequence	C-QHAASLGDLTPN
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	127343
Mouse Gene ID	140477
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	No signal obtained yet but low background observed in whole brain extracts at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA