


GOAT ANTI-DMTF1 (ISOFORM A) ANTIBODY

SKU: EB12122



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. |
| Unit Size | 100 µg |
| Storage Instructions | Aliquot and store at -20°C. Minimize freezing and thawing. |
| Synonym / Alias Names | DMTF1 cyclin D binding myb-like transcription factor 1 DMP1 DMTF hDMP1 cyclin D-binding Myb-like protein cyclin-D-binding Myb-like transcription factor 1 cyclin-D-interacting Myb-like protein 1 hDMTF1 |
| Accession ID | NP_066968.3; NP_001135798.1 |
| Blocking Peptide | EBP12122 |
| Immunogen | Peptide with sequence CPQNEADEIDSEDS, from the N-Terminus of the protein sequence according to NP_066968.3; NP_001135798.1. |
| Product Comments | This antibody is expected to recognize both reported isoform a (NP_066968.3) only. Reported variants represent identical protein: NP_066968.3, NP_001135799.1 |
| Peptide Sequence | CPQNEADEIDSEDS |
| Purification Method | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. |
| Shipping Instructions | Refrigerated |
| Predicted Species | Human, Dog, Pig, Cow |
| Reactive Species | Human, Mouse, Rat |
| Human Gene ID | 9988 |
| Mouse Gene ID | 23857 |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png |
| ELISA Detection Limit | Antibody detection limit dilution 1:16000. |
| Western Blot | Approx 70kDa band observed in Human Lung and Testis lysates (calculated MW of 74.9kDa according to NP_001135798.1). Recommended concentration: 1-3µg/ml. |
| Application Type | Pep-ELISA, WB |

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

