

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-DMTF1 (ISOFORM A) ANTIBODY

SKU: EB12122



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa


50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names DMTF1|cyclin D binding myb-like transcription factor 1|DMP1|DMTF|hDMP1|cyclin D-binding Myb-like protein|cyclin-D-binding Myb-like transcription factor 1|cyclin-D-interacting Myb-like protein 1|hDMTF1

Accession

ID

NP_066968.3; NP_001135798.1

Blocking

Peptide

EBP12122

Immunogen

Peptide with sequence CPQNEADEIDSEDS, from the N-Terminus of the protein sequence according to

NP_066968.3; NP_001135798.1.

Product

This antibody is expected to recognize both reported isoform a (NP_066968.3) only. Reported variants

Comments

represent identical protein: NP_066968.3, NP_001135799.1

Peptide

Sequence

CPQNEADEIDSEDS

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Human, Dog, Pig, Cow

Reactive

Species

Human, Mouse, Rat

9988

Human

Gene ID

Mouse

23857 **Gene ID**

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Detection

Antibody detection limit dilution 1:16000.

Limit

Western

Approx 70kDa band observed in Human Lung and Testis lysates (calculated MW of 74.9kDa according to

Blot NP_001135798.1). Recommended concentration: 1-3µg/ml.

Application

Type

Pep-ELISA, WB

Telephone: (650) 697-3600



Email: customerservice@vectorlabs.com

GALLERY IMAGES

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa