

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-DNA POLYMERASE THETA ANTIBODY

SKU: EB08503

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

DNA polymerase theta|DNA polymerase eta|PRO0327|POLH|DKFZp781A0112|polymerase (DNA directed), **Alias**

theta|POLQ **Names**

Accession

NP 955452.3

Blocking

EBP08503 **Peptide**

Peptide with sequence DTHQQLTCSNDES, from the internal region of the protein sequence according to Immunogen

NP 955452.3.

Peptide

DTHQQLTCSNDES Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human

10721 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:1000.

Limit

Blot

Western Preliminary experiments in Human Brain (Cerebellum lysates) and in lysates of cell line HeLa gave no specific

signal but low background (at antibody concentration up to 1µg/ml).

Application

Pep-ELISA **Type**