

GOAT ANTI-DOC2B ANTIBODY

SKU: EB10544

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	double C2-like domain-containing protein beta double C2 protein beta doc2-beta OTTMUSP00000006522 double C2, beta Doc2b DOC2B
Accession ID	NP_003576.2
Blocking Peptide	EBP10544
Immunogen	Peptide with sequence C-QLPVDKTEDKS, from the C-Terminus of the protein sequence according to NP_003576.2.
Peptide Sequence	C-QLPVDKTEDKS
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human, Mouse, Rat
Human Gene ID	8447
Mouse Gene ID	13447
Rat Gene ID	81820
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 48kDa band observed in Human, Mouse and Rat Heart and Kidney lysates (calculated MW of 45.9kDa according to NP_003576.2). Recommended concentration: 0.5-2µg/ml.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

