

GOAT ANTI-DOK4 ANTIBODY

SKU: EB05211

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	FLJ10488 downstream of tyrosine kinase 4 DOK4
Names	
Accession ID	NP_060580
Immunogen	Peptide with sequence C-SQGDSSEAKTPSQ, from the C Terminus of the protein sequence according to NP_060580.
Peptide Sequence	C-SQGDSSEAKTPSQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Mouse
Human Gene ID	55715
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: Approx 45-50kDa band seen in Mouse Kidney extracts. Recommended for use at 0.5-2µg/ml. Overnight incubation is recommended with this antibody.