

GOAT ANTI-DOPAMINE RECEPTOR D1 / DRD1 ANTIBODY

SKU: EB07420



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	DRD1A DADR dopamine receptor D1 DRD1
Names	
Accession ID	NP_000785.1
Blocking Peptide	EBP07420
Immunogen	Peptide with sequence C-DYDTDVSLEKIQ, from the C Terminus of the protein sequence according to NP_000785.1.
Peptide Sequence	C-DYDTDVSLEKIQ
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	1812
Mouse Gene ID	13488
Rat Gene ID	24316
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Preliminary experiments in Human Brain lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)