

GOAT ANTI-DRAK1 / STK17A ANTIBODY

SKU: EB05185



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	DAP kinase-related apoptosis-inducing protein kinase 1 DAP kinase-related apoptosis-inducing protein kinase 1 death-associated protein kinase-related 1 DRAK1 serine/threonine kinase 17a (apoptosis-inducing) STK17A
Accession ID	NP_004751.2
Immunogen	Peptide with sequence C-EEPLLQEIPGEFIY, from the C Terminus of the protein sequence according to NP_004751.2.
Peptide Sequence	C-EEPLLQEIPGEFIY
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Rabbit
Reactive Species	Human
Human Gene ID	9263
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 50kDa band seen in MOLT-4 cell lysates. Recommended for use at 1-3µg/ml. Overnight incubation is recommended with this antibody.
Application Type	Pep-ELISA, WB