

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-DRAK1 / STK17A ANTIBODY

SKU: EB05185

×

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

 $\textbf{Unit Size} \qquad 100 \; \mu \text{g}$

Storage

Instructions

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias
Names

DAP kinase-related apoptosis-inducing protein kinase 1|DAP kinase-related apoptosis-inducing protein kinase 1|death-associated protein kinase-related 1|DRAK1|serine/threonine kinase 17a (apoptosis-inducing)|STK17A

Accession

NP_004751.2

Peptide with sequence C-EEPLLQEIPGEFIY, from the C Terminus of the protein sequence according to Immunogen

NP 004751.2.

Peptide

Sequence C-EEPLLQEIPGEFIY

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions Refrigerated

Predicted

Species Human, Rabbit

Reactive

Species Human

. Human

Gene ID

Product

Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Western

Western blot: Approx 50kDa band seen in MOLT-4 cell lysates. Recommended for use at 1-3μg/ml. Overnight

Blot incubation is recommended with this antibody.

Application

Type Pep-ELISA, WB