Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-DUDULIN 2 / STEAP3 ANTIBODY

SKU: EB08378

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

tumor suppressor pHyde|six transmembrane prostate protein 3|dudulin 2|dudlin **Alias** 2|1010001D01Rik|dudlin-2|TSAP6|STMP3|STEAP family member 3|STEAP3 **Names**

Accession

NP_878919.2; NP_001008410.1; NP_060704.2

Blocking

ID

EBP08378 **Peptide**

Peptide with sequence C-QPYVQESQNKFFK, from the internal region of the protein sequence according to **Immunogen**

NP_878919.2; NP_001008410.1; NP_060704.2.

Product This antibody is expected to recognise all reported isoforms (NP_878919.2; NP_001008410.1 and

Comments NP_060704.2).

Peptide

C-QPYVQESQNKFFK Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Instructions

Refrigerated

Predicted

Human **Species**

Reactive

Human

Species

Human 55240

Gene ID **Product**

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Antibody detection limit dilution 1:16000.

Limit

Blot

Western Blot: Approx 50-55kDa band observed in Human Liver and Pancreas lysates (calculated MW of Western 55.7kDa according to NP_878919.2 and of 54.6kDa according to NP_001008410.1 and NP_060704.2).

Recommended concentration: 1-3µg/ml.





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Application Pep-ELISA, WB