


GOAT ANTI-DUSP14 / MKP6 ANTIBODY

SKU: EB05320



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	DUSP14 dual specificity phosphatase 14 MKP6 MKP-L MKP-1 like protein tyrosine phosphatase OTTHUMP00000164064 OTTHUMP00000164065
Accession ID	NP_008957.1
Blocking Peptide	EBP05320
Immunogen	Peptide with sequence CEKESRHLMPYWGI, from the C Terminus of the protein sequence according to NP_008957.1.
Peptide Sequence	CEKESRHLMPYWGI
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated

Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	11072
Mouse Gene ID	56405
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Testis. Recommended concentration: 3.5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	No signal obtained yet but low background observed in Human Liver extracts at up to 1µg/ml.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

