

GOAT ANTI-DUSP8 ANTIBODY

SKU: EB05235



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	serine/threonine specific protein phosphatase H1 phosphatase, vaccinia virus homolog HVVH-5 HVVH8 HB5 dual specificity phosphatase 8 DUSP8
Accession ID	NP_004411
Immunogen	Peptide with sequence *AGDRLPRKVMDAK-C, from the N Terminus of the protein sequence according to NP_004411.
Peptide Sequence	*AGDRLPRKVMDAK-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	1850
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Western blot: No signal obtained yet but low background observed in Human Brain extracts at upto 2µg/ml.