



Telephone: (650) 697-3600



GOAT ANTI-DUSP8 ANTIBODY

SKU: EB05235

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / serine/threonine specific protein phosphatase|H1 phosphatase, vaccinia virus homolog|HVH-5|HVH8|HB5|dual

Alias specificity phosphatase 8|DUSP8

Names

Accession NP 004411

Peptide with sequence *AGDRLPRKVMDAK-C, from the N Terminus of the protein sequence according to Immunogen NP_004411.

Peptide

*AGDRLPRKVMDAK-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse **Species**

Human

1850 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Antibody detection limit dilution 1:1000. **Detection**

Limit

Western Blot

Western blot: No signal obtained yet but low background observed in Human Brain extracts at upto 2μg/ml.